

Greencare[®]

by  ortho
europe[®]



Ep1c

Fully-adaptable wheelbase

- Double brace frame fashioned from aircraft grade aluminum.
- Designed primarily for indoor use.
- Dynamic and lightweight.

ortho-europe.com
greencaremobility.com

Greencare[®]

by **ortho**
europe[®]

Epic

This lightweight, powered wheelchair uses the highly acclaimed Greencare standard DBI double brace frame fashioned from aircraft grade aluminium, which is strong yet light. It is primarily designed for indoor use but also has capability for some outdoor use.

Standard Features

- Greentyre puncture proof tyres
- Ultra soft easy-to-clean seat & backrest upholstery
- 10" / 250mm high traction rear wheels
- 8" / 200mm large profile front wheels
- Detachable batteries
- Linx Dynamic control system
- Joystick position variation
- Attendant or dual control options
- Greencare lap belt
- 10° folding / non-folding back posts
- Safety horn
- LED colour battery Indicator (on joystick)
- Double cross brace
- Memory foam cushion

Features

Greentyre
puncture proof tyres

Linx Dynamic
control system

Approved
for vehicle transportation



Specification

Description	Measurement
Seat width	15", 17", 18" / 380mm, 430mm, 460mm
Seat depth	15" or 17" / 380mm or 430mm
Push handle height	35.5" / 900mm
Backrest height	19" / 480mm
Seat to ground height	Standard
Weight capacity	Rated for users up to 125kg
Maximum speed	4 mph / 6.4 km/h
Maximum safe slope	8°
Motor power	200W
Weight	28kg (without batteries) / 41kg (with batteries)
Turning circle	1640mm - 1800mm
Battery life cycle unladen	17.4 miles / 28km
Battery life cycle under load	7.7 miles / 12.5km
Warranty	24 months against manufacturing defects



Ortho Europe

Ability House, Nuffield Way,
Abingdon, OX14 1RL
+44 (0)1235 552 895
info@ortho-europe.com

Scan the QR code
to find out more about
the Epic:

