

# Greencare<sup>®</sup>

by  ortho  
europe<sup>®</sup>



## DB1 Attendant

Fully-adaptable wheelbase

- Configurable chair to meet user needs
- Lightweight design
- Easily maneuverable with a choice of rear wheels for stability and performance

[ortho-europe.com](http://ortho-europe.com)  
[greencaremobility.com](http://greencaremobility.com)



by  ortho  
europe®

## DB1 Attendant

Made from an aluminium double brace for added strength and rated for users up to 125kg, the DB1 Attendant is designed to excel in rear pushing ability.

The purpose of the DB1 Attendant is to enable the occupant to relax whilst being pushed at the rear of the chair. This high comfort, high quality wheelchair is reliable, durable, strong and can be modified in any way to suit the user's exact needs. It is one of the lightest wheelchairs on the market with this feature set and can be easily folded and stowed for transportation in a small car.

The DB1 Attendant satisfies the requirements of ISO 7616-19.

### Standard specification

- Greentyre puncture proof tyres
- Ultra soft easy-to-clean seat & backrest upholstery
- Flip-back armrests
- Folding backrest
- Removable swing-away legrests
- Double cross brace
- Plastic side panels
- Memory foam cushion
- 110° swing-away footrests
- 10° backrest
- 12" or 16" rear wheels
- 7.5" / 190mm front castor wheels



### Features

**Greentyre**  
puncture proof tyres

**Ultra soft easy-to-clean**  
seat & backrest upholstery

**High comfort, high quality**  
reliable, durable and strong

# Specification

Description	Measurement
Seat width	15" – 20" / 380mm – 508mm
Seat depth	115", 17" or 19" / 380mm, 430mm or 480mm
Standard push handle height	35.5" / 900mm
Seat back height	19" / 480mm
Seat to ground height	19" / 480mm
Weight capacity	Rated for users up to 125kg
Warranty	24 months against manufacturing defects
Crash testing	ISO 7616-19





**Ortho Europe**

Ability House, Nuffield Way,  
Abingdon, OX14 1RL  
+44 (0)1235 552 895  
[info@ortho-europe.com](mailto:info@ortho-europe.com)

**Scan the QR code**  
to find out more about  
the Attendant DBI:

