

Greencare[®]

b:



DB1 Junior

Fully-adaptable wheelbase

- Made from both a single or double aluminium brace for choice in weight and added strength
- Easily-adaptable to accommodate growth or changing needs
- Coloured Greentyre puncture proof tyres

ortho-europe.com
greencaremobility.com



by **ortho**
europe®

DB1 Junior

The DB1 Junior wheelchair is an extremely lightweight, simple to use, and easy to transport wheelchair that has been designed to accommodate changing user needs.

Designed to move easily within its environment, the DB1 Junior features multi adjustable wheels that can be easily moved to ensure the user maximises all their strength to get the most out of the chair. The DB1 Junior can also easily be adapted to accommodate changing needs as the user grows.

Standard Features

- Coloured Greentyre puncture proof tyres
- Ultra soft, easy-to-clean seat & backrest upholstery
- Flip-back armrests
- Folding backrest
- Removable leg rests
- Single cross brace
- Memory foam cushion
- Coloured polyurethane hand rims
- Custom coloured spoke guards
- 90° swing-away footrests
- 3 / 10° backrest
- 20", 22" , 24" quick release self propelling wheels
- 6" / 150mm castor wheels

Features

Maintenance-free
puncture proof tyres

Ultra soft easy-to-clean
seat & backrest upholstery

Easily refurbished
to grow with the user



Specification

Description	Measurement
Seat width	11" - 15" / 280mm, 380mm
Seat depth	13" or 15" / 330mm or 380mm
Standard push handle height	35.5" / 900mm
Seat back height	19" / 480mm
Seat to ground height	17.5" / 450mm
Weight capacity	Rated for users up to 100kg
Warranty	24 months against manufacturing defects





Ortho Europe

Ability House, Nuffield Way,
Abingdon, OX14 1RL
+44 (0)1235 552 895
info@ortho-europe.com

Scan the QR code
to find out more about
the DBI Junior:

