



by  **ortho**  
europe®



# Pelvic positioning belts

For postural support and stability

- Designed to provide postural support and stability to the pelvis of a wheelchair user
- A range of two- and four-point belts

[ortho-europe.com](http://ortho-europe.com)  
[greencaremobility.com](http://greencaremobility.com)



by  **ortho**  
europe®

## Pelvic positioning belts

The SeatMatch pelvic belt range is designed to provide postural support and stability to the pelvis of a wheelchair user. Packs of mounting components are available, ensuring no compromise when fitting any of our pelvic belts.

These belts must not be used as the sole method of occupant restraint during vehicle transportation.

We provide a range of two- and four-point belts, both with either a single- or dual-pull capability. Our metal, single-pull buckle is our heavy-duty option, whilst our slimline plastic buckle has dual-pull side-release buttons.

### Features

**Provides postural support** and stability to the pelvis

**Two- and four-point belts** both with either a single- or dual-pull capability

**Slimline plastic buckles** with dual-pull, side-release buttons



## Specification - metal buckles

Two-point single pull - padded



Four-point single pull - padded



Metal buckle	Small	Medium	Large
Belt size	150mm x 44mm	150mm x 60mm	150mm x 70mm
Webbing	25mm	38mm	50mm
Two-point single pull - padded	HD2001	HD2002	HD2003
Four-point single pull - padded	HD5002	HD5003	HD5004

## Specification - plastic buckles

Two-point dual pull - padded



Four-point dual pull - padded



Plastic buckle	Small	Medium	Large
Belt size	150mm x 44mm	150mm x 60mm	150mm x 70mm
Webbing	25mm	38mm	50mm
Two-point single pull - padded	HD2002S	HD2003S	HD2004S
Two-point dual pull - padded	HD3002S	HD3003S	HD3004S
Four-point dual pull - padded	HD5002S	HD5003S	HD5004S

## Pipe clamp (for 19, 22 and 25mm tubes)

Part number	Description
E12 023	Pipe clamp set with 2 clamps and 4 reduction rings





**Ortho Europe**

Ability House, Nuffield Way,  
Abingdon, OX14 1RL  
+44 (0)1235 552 895  
info@ortho-europe.com

**Scan the QR code**  
to find out more about  
the pelvic positioning  
belts:

