



by  **ortho**
europe®



Harnesses

- A range of harnesses to provide support and stability to the upper trunk
- Designed to be used with a pelvic belt

ortho-europe.com
greencaremobility.com



by  **ortho**
europe®

Harnesses

Designed to provide support of the torso for improved overall function and positioning, Seatmatch postural supports are available in three support styles and with either structured or dynamic construction to accommodate clients with varying positioning needs and levels of trunk control.

The main characteristics are:

- Adjustable-length top straps, with plastic pull-tight buckles and cam-lock features
- Lower straps are fixed-length, with angle-adjusting plastic buckles
- Complete with Metal triple-slot mounting plates
- Nylon-covered neoprene pad can be easily removed for machine washing
- Each type is available in three sizes

Features

Adjustable-length straps

Nylon-covered neoprene pad
easily removed for machine washing



Lower straps
with angle-adjusting plastic buckles



Specification


Dynamic Harness

| Part No | HU3002 | HU3003 | HU3004 |
|---------|--------|--------|--------|
| | Small | Medium | Large |
| A | 413mm | 457mm | 514mm |
| B | 241mm | 267mm | 298mm |
| C | 246mm | 280mm | 317mm |
| D | 50mm | 50mm | 50mm |



Sternum Harness

| Part No | HU4002 | HU4003 | HU4004 |
|---------|--------|--------|--------|
| | Small | Medium | Large |
| A | 394mm | 470mm | 527mm |
| B | 235mm | 235mm | 235mm |
| C | 235mm | 235mm | 235mm |
| D | 95mm | 114mm | 127mm |
| E | 374mm | 444mm | 470mm |



Shoulder Harness

| Part No | HU5002 | HU5003 | HU5004 |
|---------|--------|--------|--------|
| | Small | Medium | Large |
| A | 254mm | 280mm | 317mm |
| B | 89mm | 102mm | 102mm |
| C | 44mm | 51mm | 51mm |





Ortho Europe

Ability House, Nuffield Way,
Abingdon, OX14 1RL
+44 (0)1235 552 895
info@ortho-europe.com

Scan the QR code
to find out more about
the harnesses:

